Laboratoriemedicin

Godkänt datum 2026-03-31

P-Faktor IX (koag)

Bakgrund

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Nivå /CV%: 0,1/20

Referenslitteratur

1. Laurells Klinisk kemi i praktisk medicin Studentlitteratur 2018, Koagulationsrubbnings p.171-207.
2. 6th Edition of the European Pharmacopoeia, General Chapter 5,3 Statistical analysis of results of biological assays and tests.
3. Gray E, Tubbs J, Thomas S et al. Measurement of activated Factor IX in Factor IX concentrates: Correlation with I Vivo thrombogenicity. *Thromb Haemost* 1995; 73 (4): 657-9
4. Pickering W M, Gray E. The effect of activated Factor IX on the Factor IX coagulant and NAPTT activity of high-purity Factor IX concentrates. *J Thromb Haemost* 2007; 5, Supplement 2: P-T-156.
5. Hubbard AR, Dodt J, Lee T, et al. Recommendations on the potency labelling of factor VIII and factor FIX concentrates. *J Thromb Haemost* 2013; 11:988-9.
6. Verifiering Sysmex CS-5100, Specialkoagulation, Malmö (DokumentID 32478042)
7. Certifikat till Factor IX Deficient Plasma, (George King Biomedic Inc, USA)
8. Bipacksedel till Dade® Actin® FS Aktiverad PTT-reagens. 11528763_sv, Senaste versionen (Siemens Healthineers).
9. Sysmex CS-5100 System Reference Guide, Senaste versionen (Siemens Healthineers).
10. CS-5100 System Evaluation and Check Algorithm (OUS), Senaste versionen (Sysmex).