Parasiter Malaria

Bakgrund

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Referenslitteratur

- WHO. World Malaria Report 2021. Geneva: World Health Organisation; 2021.
- WHO. WHO Guidelines for Malaria. Geneva: World Health Organisation; 2022.
- WHO. Malaria rapid diagnostic test performance: results of WHO product testing of malaria RDTs: round 8 (2016–2018). Geneva: World Health Organisation; 2018.
- WHO. WHO Evidence Review Group on Malaria Diagnosis in Low Transmission Settings Meeting report. Geneva: WHO; 2013.
- Morris U, Aydin-Schmidt B. Performance and Application of Commercially Available Loop-Mediated Isothermal Amplification (LAMP) Kits in Malaria Endemic and Non-Endemic Settings. *Diagnostics*. 2021;11(2):336.
- Oyegoke OO, Maharaj L, Akoniyan OP, Kwoji I, Roux AT, Adewumi TS, et al. Malaria diagnostic methods with the elimination goal in view. *Parasitology Research*. 2022.
- Ljolje D, Abdallah R, Lucchi NW. Detection of malaria parasites in samples from returning US travelers using the Alethia® Malaria Plus LAMP assay. *BMC Research Notes*. 2021;14(1).
- Selvarajah D, Naing C, Htet NH, Mak JW. Loop-mediated isothermal amplification (LAMP) test for diagnosis of uncomplicated malaria in

endemic areas: a meta-analysis of diagnostic test accuracy. *Malaria Journal*. 2020;19(1).

- Rypien C, Chow B, Chan WW, Church DL, Pillai DR. Detection of Plasmodium Infection by the illumigene Malaria Assay Compared to Reference Microscopy and Real-Time PCR. *Journal of Clinical Microbiology*. 2017;55(10):3037-45.
- Hartmeyer GN, Hoegh SV, Skov MN, Kemp M. Use of Loop-Mediated Isothermal Amplification in a Resource-Saving Strategy for Primary Malaria Screening in a Non-Endemic Setting. *The American Journal of Tropical Medicine and Hygiene*. 2019;100(3):566-71.
- Polley SD, González IJ, Mohamed D, Daly R, Bowers K, Watson J, et al. Clinical Evaluation of a Loop-Mediated Amplification Kit for Diagnosis of Imported Malaria. *The Journal of Infectious Diseases*. 2013;208(4):637-44.
- Cuadros J, Martin Ramírez A, González IJ, Ding XC, Perez Tanoira R, Rojo-Marcos G, et al. LAMP kit for diagnosis of non-falciparum malaria in Plasmodium ovale infected patients. *Malaria Journal*. 2017;16(1).
- Charpentier E, Benichou E, Pagès A, Chauvin P, Fillaux J, Valentin A, et al. Performance evaluation of different strategies based on microscopy techniques, rapid diagnostic test and molecular loop-mediated isothermal amplification assay for the diagnosis of imported malaria. *Clinical Microbiology and Infection*. 2020;26(1):115-21.
- Nolasco O, Infante B, Contreras-Mancilla J, Incardona S, Ding XC, Gamboa D, et al. Diagnosis of Plasmodium vivax by Loop-Mediated Isothermal Amplification in Febrile Patient Samples from Loreto, Perú. *The American Journal of Tropical Medicine and Hygiene*. 2020;103(4):1549-52.
- Piera KA, Aziz A, William T, Bell D, González IJ, Barber BE, et al. Detection of Plasmodium knowlesi, Plasmodium falciparum and Plasmodium vivax using loop-mediated isothermal amplification (LAMP) in a co-endemic area in Malaysia. *Malaria Journal*. 2017;16(1).
- Meridian Bioscience. Alethia Malaria and Malaria PLUS DNA Amplification Assays, User Information. Cincinnati, Ohio, USA: Meridian Bioscience; 2020.
- Egen validering/verifiering av metoden genomfördes 2021-2022.