

Instructions for submitting DNA sequencing samples

Sample labelling

If you submit your samples in tubes: Use your initials and a sample number, e.g. PS01, PS02.

Sample names should not contain any special characters (! # & ?, . space or letters Å, Ä, Ö).

Mark your tubes clearly with the sample name.

If you have samples in 96-plates: Mark the plate with your initials or name.

Sample containers

Dried samples should be delivered in the following types of containers:

96 well ABI-PCR plate, correctly sealed

200 µL PCR tube

1.5 mL standard Eppendorf tube

0.5 mL standard Eppendorf tube

You can leave in your samples

Monday to Friday: 08.30 am to 04.00 pm

The ABI 3730 will be run on Tuesdays and Thursdays. The deadline for same day processing is 08.30 am.

Delivery address

Please fill in the DNA-sequencing order form and send it together with your samples to the following address:

RSKC, Wallenberglaboratoriet vid Universitetssjukhuset MAS,
ingång 46, plan 3, 205 02 Malmö

Please send also a copy of the order form to:

dnasekvens.rskc.umas@skane.se

tfn: 040-33 66 34 / 33 10 53

Supported chemistry

We can analyze reactions run by the Applied Biosystems BigDyeTerminator v 3.1 chemistry.

Reporting of results

Your results will be delivered via E-mail to the address specified in contact details.

Price

See separate pricelist.

DNA-Sequencing Order form

Date/Datum	
Name/Namn	
Department/Avdelning	
Cost Center/Kostnadställe	
Invoice/Faktura address	
E-mail	
Phone/Tfn	

Tube ID/Well	Sample name
A01	
B01	
C01	
D01	
E01	
F01	
G01	
H01	
A02	
B02	
C02	
D02	
E02	
F02	
G02	
H02	
A03	
B03	
C03	
D03	
E03	
F03	
G03	
H03	
A04	
B04	
C04	
D04	
E04	
F04	
G04	
H04	

A05	
B05	
C05	
D05	
E05	
F05	
G05	
H05	
A06	
B06	
C06	
D06	
E06	
F06	
G06	
H06	
A07	
B07	
C07	
D07	
E07	
F07	
G07	
H07	
A08	
B08	
C08	
D08	
E08	
F08	
G08	
H08	
A09	
B09	
C09	
D09	
E09	
F09	
G09	
H09	
A10	
B10	

C10	
D10	
E10	
F10	
G10	
H10	
A11	
B11	
C11	
D11	
E11	
F11	
G11	
H11	
A12	
B12	
C12	
D12	
E12	
F12	
G12	
H12	